

# Dialog+ Accessories

Convenient. Flexible. Efficient.



Hemodialysis

# Dialog+ Accessories

Sophisticated extensions for an efficient workflow






**1 Universal Front Tray\***  
Protects front components against any kind of damage and serves as a storage area for syringes or other small parts.



**5 Sample Port**  
The dialysis fluid sample port is integrated into the dialyzer inlet tubing. Direct, needle-free sampling of dialysis fluid is possible with a syringe.



**2 Monitor Mini-Shelf\***  
Removable Monitor Mini-Shelf for small parts such as syringes or clamps. Located below the monitor.



**6 Crank/Crank Pin**  
Use the crank or crank pin to manually operate the blood pump in case of power failure.



**3 Dialyzer Holder, Metal**  
Swivel rotatable aluminum dialyzer holder, including holder for bloodlines.



**7 Combi Shelf Holder**  
Removable metal tray that creates a flat surface behind the monitor. Helps to place box of disposables and other devices like Perfusor®, Infusomat®, etc. Features a built-in recessed holder for small parts.



**4 Rinse Bucket**  
Removable rinse bucket to collect priming saline liquid, with tubing hold.



**8 Multifunctional Tray**  
Plastic tray that enables flexible storage of any kind of disposables, devices, and other objects. Easy to remove and clean anytime, without additional tools.

\* **1 2** Not available in combination with each other.



**9 Box Comfort**  
Multifunctional box with a slot for storing patient records and a separate basket to hold all the disposables.



**13 Disinfection Canister Holder**  
Stand mounted on the rear panel of the Dialog+ machine holds canisters containing 6 or 10 liters of disinfecting agent. Can be easily removed for cleaning.



**10 Holder or Basket for Blood Pressure Cuff\*\***  
Hang the blood pressure cuff on this convenient holder on the side of the dialysis machine. Store the blood pressure cuff in a basket on left side of the dialysis machine.



**14 Universal Storage Tray**  
Tray on rear panel for storing a single-site osmosis system or, alternatively, two 10-liter canisters of disinfecting agent. Only available for Dialog+ standard machines.



**11 Concentrate-Bag Holder\*\***  
Convenient holder for hemodialysis concentrates in bags.



**12 Record Holder\*\***  
Keeps the patient record and documents beside the machine for easy access.

\*\* 10 11 12 Not available in combination with each other.

# Product Specifications

	Product Code	Product Description	Packaging Unit
1	7105239	Universal Front Tray	1 unit
2	7102872	Monitor Mini-Shelf	1 unit
3	7107426	Dialyzer Holder, Metal	1 unit
4	7105237	Rinse Bucket	1 unit
5	7102876	Sample Port, AP Version	1 unit
5	7102867	Sample Port, ROW Version	1 unit
6	34560157	Crank Pin for the Blood Pump	1 unit
6	34518355	Crank for the Blood Pump	1 unit
7	7102871	Combi Shelf Holder (SW Versions < 9.xx)	1 unit
7	7102890	Combi Shelf Holder (SW Versions ≥ 9.xx)	1 unit
7	7102869	Non-Skid Padding for the Combi Shelf Holder	1 unit
8	7105238	Multifunctional Tray	1 unit
9	7107322	Box Comfort	3 units plus 4 red, 4 blue, and 4 lucent rings
10	7102781	Holder for Blood Pressure Cuff	1 unit
10	7102865	Basket for Blood Pressure Cuff	1 unit
10	7102372	ABPM Cuff, Latex-Free, incl. Tubing, Small	1 unit
10	7102771	ABPM Cuff, Latex-Free, incl. Tubing, Medium	1 unit
10	7102380	ABPM Cuff, Latex-Free, incl. Tubing, Large	1 unit
11	7102877	Concentrate-Bag Holder	1 unit
12	7102873	Record Holder	1 unit
13	7102277	Disinfection Canister Holder (for 6 l or 10 l)	1 unit
14	7105500	Universal Storage Tray	1 unit

