

# Comfort-4

# Plus

# THERAPY CHAIR

HIGHEST COMFORT, ATTRACTIVE DESIGN IN ADDITION TO FUNCTIONALITY AND SAFETY

- The Comfort-4 Therapy Chairs are designed to assist and facilitate the work of physicians and medical staff
- They increase the patient's comfort considerably and therefore help to relieve tension during treatment
- Intelligent Control System with microprocessor
- Several optional accessories for optimal choice



PROFESSIONAL SEATING SOLUTIONS BY

**Digiterm**<sup>®</sup>  
medical equipment

# COMFORT-4 PLUS THERAPY CHAIR

ART.-NR. 7130368

Optional accessories	Article Nr.
Comfort-T1 Bed-side Table	6200501
Transparent protecting plastic cover for the legrest	7130240
4 channel audio unit (built in the control box)	7130400
Central operating brake with 4x diameter 10 cm twin wheel castors	7130402
Motorised footrest	7130404
Headrest pillow with built-in stereo speakers (only together with 4 channel audio unit, Art.-Nr. 7130400)	7130405
Reading lamp with flexible arm	7130408
SSP motor control for smooth operation (Soft Start-Stop Process)	7130409
Central operating total and directional lock with 4x diameter 12,5 cm twin wheel castors	7130413
Central operating total and directional lock with 4x diameter: 15 cm twin wheel castors	7130414
Stainless steel armrest protecting frame (one pair)	7130440
Plastic side cover on both sides of seat (handles on both sides included)	7130450
Trendelenburg foot-switches on both sides of underframe (no underframe plastic cover included)	7130482
Terry cloth (Frottier) protecting cover for upholstered parts (standard white colour)	7130501
Infusion bottle holder (IV pole) with height adjustment	7130506
Washable net bag for storing Patient's belongings	7130507
Paper roll holder	7130532
Integrated multifunctional table (built on the chair)	7130600
Adjustable monitor arm for mini TV or DVD player (mounting kit with antennae and power cables hidden in the arm support) (only together with Built-in Medical Power supply, Art.-Nr. 7130611)	7130610
Built-in Medical Power Supply for LCD mini TV or DVD player	7130611
Rechargeable battery back-up operation in case power failure	7130650
Individual colour for Standard shape upholstery, selectable from more than 40 diff. colours SKAI-TUNDRA	7130690
Special flame retardant high-tech artificial leather for Standard shape upholstery, selectable from more than 15 diff. colours PALMA NF (meets M1 and some other Standards) no additional colour surcharge	7130692
3-Dimensional shape upholstery design in combination with special flame retardant high-tech artificial leather, selectable from more than 15 diff. colours PALMA NF (meets M1 and some other Standards) no additional colour surcharge	7130694
Distributor's individual logo on top of backrest	7130695
3-Dimensional shape upholstery design that makes Standard shape upholstery even more comfortable and attractive. Available in single or dual-colour combination, selectable from more than 40 diff. colours TUNDRA no additional colour surcharge	7130700
Flexible height adjustable armrest for universal medical applications (one pair)	7130800
Integrated gas spring in the arm support for smoother and safer operation (one pair)	7130810
Push bar (stainless steel handle)	7130811
Removable footrest for standalone exercise applications	7130812

## General technical data

Seat height adjustment	54 - 76 cm	Size of castors	4x diameter 10 cm swivel castors with total lock
Backrest adjusting range	80° (-13° under horizontal in Trendelenburg position)	Weight of chair	80 kg
Seat tilt adjustment	20°	Maximal load	200 kg weight capacity
Legrest adjustment	30°	Upholstery colour	Standard light-grey (PERLE) and more than 40 different colours upon request
Footrest adjustment	29 cm	Standard Power supply	230 V, 50/60 Hz (110 V, 127 V optional)
Dimensions in seat position (L x W)	145 x 92 or 145 x 80 cm *	Motor's Operation Voltage	24V DC, IP66
Dimensions in full recline position (L x W)	min. 215 x 80, max. 227 x 92 cm *	Degree of protection for motors	Protected against splashing water - IP66
Size of upholstery	L 215 x W 60 x H 9.5 cm	Applied directives and standards	IEC 601-1 (EN60601-1 in Europe), IEC 601-2-52 (EN60601-2-52 in Europe); 93/42/EEC; CE

\*Depend on armrest and footrest setting

EN. 12/2 Nr.7081600



**Digiterm Ltd.**  
 Ipar út. 7., H-3200 Gyöngyös  
 Tel.: +36 37 311 990  
 Fax: +36 37 311 893  
 E-mail: sales@digiterm.hu  
 www.digiterm.hu



Approved by TÜV Product Safety

