


HIGHEST COMFORT, ATTRACTIVE DESIGN IN ADDITION TO FUNCTIONALITY AND SAFETY

- The Comfort-T1 Bed-side Tables are designed to assist and facilitate the work of physicians and medical staff
- They increase the patient's comfort considerably during long-term treatments
- Practical complementary equipment for all type of Comfort Chairs



# COMFORT-T1 BED-SIDE TABLE

ART.-NR. 6200501



\* Comfort-3 Classic Therapy Chair with Comfort-T1 Bed-side table

\* Comfort-3 Classic Therapy Chair with Comfort-T1 Eco Bed-side table

## Comfort-T1 Bed-side table

Suitable for any type of Comfort Therapy Chair and can usefully complement the latter during treatment.

It can be used on both sides of the chair. The height of the table can be adjusted as required using the lever on its support stand.

The trays can be easily cleaned and lifted out. The table can be tilted at an appropriate angle.

## Comfort-T1 Eco Bed-side table

It is economical version without the tilting function. Article Nr. 6200501E

## General technical data

Height adjustment	75-98 cm
Size of trays	40 x 56 cm and 40 x 17 cm
Table tilt angle	0±15° ±30°
Frame colour	RAL 9001
Weight	15 kg
Maximum load	15 kg



### Digiterm Ltd.

Egri u. 30., H-3200 Gyongyos

Tel.: +36 37 311 990

Fax: +36 37 311 893

E-mail: sales@digiterm.hu

www.digiterm.hu

Approved by TÜV Product Service

